

COVID-19/FLU A & B POINT OF CARE ANTIGEN TEST

Point-of-Care (POC), Status™ Covid-19/Flu A & B allows effective screening of COVID-19 infection on a large scale.



Product Features:

- Market leading accuracy
 - COVID-19 - Sensitivity 93.9%, Specificity 100%
 - Flu A - Sensitivity 91.4%, Specificity 95.7%
 - Flu B - Sensitivity 87.6%, Specificity 95.9%
- Visually read in 15 minutes
- Flocked nasopharyngeal swab for superior specimen collection and patient comfort
- FDA Emergency Use Authorization (EUA)



Order Information:

ORDERING INFORMATION		
PRODUCT	# OF TESTS	ITEM #
Status™ Covid-19/Flu A&B	25	33225

This product has not been FDA cleared or approved, but has been authorized by FDA under an Emergency Use Authorization (EUA) for use by laboratories certified under CLIA that meet the requirements to perform moderate, high or waived complexity tests. This test is authorized for use at the Point-of-Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation. This product has been authorized only for the detection of proteins from SARS-CoV-2 influenza A and influenza B, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

CORREA LOGISTICS

YOUR SINGLE SOURCE TEST KITS PARTNER

CORREA LOGISTICS OFFERS:

- ✓ Unmatched Range of Covid-19 Testing Options
- ✓ Shipping within 24 hours

OTHER COVID-19 TESTING OPTIONS




CareStart™ COVID-19 RAPID ANTIGEN TEST

- Intended Point-of-Care (POC) designated test with a 10-minute processing time, CareStart™ COVID-19 Antigen Test allows effective screening of COVID-19 infection on a large scale.

FaStep™ COVID-19 RAPID TEST SEROLOGY KIT

- Detection of IgM indicated recent infection and early diagnosis
- First CLIA-Waived test with an EUA
- 10-minute Test



Healgen™ COVID-19 RAPID TEST SEROLOGY KIT

- Multi-variable analysis of immunoglobulin IgG & IgM
- Emergency Use Authority (EUA) by the FDA
- 10-minute test

MD-Bio™ COVID-19 RAPID MOLECULAR TEST

- The Bio BCC19 Test Kit is a reverse transcriptase loop-mediated isothermal amplification (RT-LAMP) test
- Operating at the speed of antigen testing generating results in 30 minutes
- Up to 96 samples can be run at a time with one unit



Disclaimers:

- This test has been authorized by FDA under an EUA for use by authorized laboratories.
- This product is intended for professional use (POC) by medical professionals
- Not for the screening of donated blood.